



QUANSYS


B I O S C I E N C E S



Competitive and Sandwich ELISAs

Sandwich ELISAs

- ▶ Primary or capture antibody is applied to the well.
- ▶ The sample containing the antibody's target (antigen) is added.
- ▶ The complimentary or match pair antibody which is labeled is added to the well.
- ▶ If the antigen was present in the sample, it has then bound to the capture antibody and hence the detection binds, building the antibody/antigen complex.
- ▶ The labeled detection will then be read or detected.




Sandwich ELISAs

▶ Pros:

- ▶ Sensitive
- ▶ Standardized method


▶ Cons:

- ▶ More steps in Assay
- ▶ Two antibodies
- ▶ Only applicable on large molecules



Competitive ELISAs

- ▶ Primary or capture antibody is applied to the well.
- ▶ The sample containing the antibody's target (antigen) is added.
- ▶ In the next step, labeled antigen is added.
- ▶ Whatever binding sites remain, they will be occupied by the labeled antigen
- ▶ The detector (SHRP) is added.
- ▶ Where no antigen was present in the sample, the labeled antigen will have bound and create a signal.
- ▶ If the sample contained antigen, no labeled antigen will bind, causing no signal. -inverse curve




Competitive ELISAs

▶ Pros:

- ▶ Needs one antibody
- ▶ Small molecules

▶ Cons:

- ▶ Less sensitive



Development Strategy

- ▶ Solution: Combine both sandwich and competitive ELISA in the same panel of assays so a single dilution is used.
- ▶ Benefits:
 - ▶ All assays are within one well.
 - ▶ Panels with different ranges (ug/ml and ng/ml) can be multiplexed by converting sandwich ELISAs to competitive to alleviate multiple dilutions
 - ▶ When only one antibody is available for an assay in a panel, the assay no longer becomes a show-stopper
 - ▶ Small and large molecule assays can still be in the same panel

CASE STUDY # 1

SFU


- ▶ Validation of Custom Multiplex Array Against Individual ELISAs
 - ▶ Simon Fraser University, BC, Canada
 - ▶ Female Reproductive Health: Urine samples from Guatemala
 - ▶ Adiponectin, Cortisol, E1G, FSHb, HCGb, and C-Peptide
 - ▶ 4 sandwich ELISAs and 2 competitive ELISAs in one well
 - ▶ Testing in parallel to Bayer ACS:180 Clinical Analyzer
 - ▶ Pearson Correlation Coefficient: (≥ 0.75)
 - ▶ [Am J Hum Biol.](#) 2012 Jan-Feb;24(1):81-6. doi: 002/ajhb.21229. Epub 2011 Nov 28.



CASE STUDY # 1



► Correlation between two methodologies



CASE STUDY # 1



▶ Assay Performance: Sensitivity and Reproducibility

	Quansys Multiplex			Traditional ELISAs		
	Sensitivity	Intra-Assay CV	Inter-Assay CV	Sensitivity	Intra-Assay CV	Inter-Assay CV
Adiponectin	0.023 ng/ml	10%	6.90%	0.156 ng/ml	4.4%	6.2%
Free Cortisol	0.343 ng/ml	7.30%	8.50%	0.057 ng/ml	10.5%	13.4%
C-Peptide	0.090 ng/ml	9.30%	7.70%	2 ng/ml	3.9%	8.5%
E1G	0.252 ng/ml	9.70%	8.20%	1.45 ng/ml	7.9%	8.5%
FSHb	0.017 ng/ml	7.20%	7.30%	0.143 ng/ml	3.8%	6.5%
HCGb	0.035 ng/ml	7.10%	7.50%	0.003 ng/ml	3.5%	5.8%

- ▶ Summary: “This multiplex technology provides a more economic, rapid, and ecologically sound alternative to individual assays for studies requiring the measurement of multiple biomarkers per biospecimen.”

Questions?

Thank You!